



Unilong Industry Co., Ltd

Add: No.2000 Shunhua Rd, High-Tech Zone , Jinan City, Shandong Province, China


Contact: +86-18653132120

Website: www.unilongmaterial.com

Product name: Alizarin Red S

CAS: 130-22-3

Molecular weight: 342.26



Item	Standard
Strength	100%
coloured light	Approximate to micro
Moisture content (%)	≤5%
Water insoluble matter (%)	≤0.5%
Fineness (um)	≤5%



Unilong Industry Co., Ltd

Add: No.2000 Shunhua Rd, High-Tech Zone , Jinan City, Shandong Province, China

Contact: +86-18653132120

Website:www.unilongmaterial.com

Application:

Alizarin red S can form colored compounds with many metal ions, which can be used for color reaction and colorimetric determination of zirconium, thorium, aluminum, titanium and beryllium. It can be used as microscopic staining agent, in vivo staining of nerve tissue, chromosomal staining in plant cytology and reagents for the determination of belladonna base, as well as for color matching in woolen, worsted, carpets and blankets.

Package:

25kg/drum

Unilong Industry